

RHBDF1 Polyclonal Antibody

Catalog No: #30163



Package Size: #30163-1 50ul #30163-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RHBDF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human RHBDF1 (NP_071895.3).
Other Names	C16orf8, Dist1, EGFR-RS, gene-89, gene-90, hDist1
Accession No.	Swiss-Prot#:Q96CC6NCBI Gene ID:64285
Uniprot	Q96CC6
GeneID	64285;
Calculated MW	97kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using RHBDF1 at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using RHBDF1 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only