

FNDC4 Polyclonal Antibody

Catalog No: #30164



Package Size: #30164-1 50ul #30164-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FNDC4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IHC
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human FNDC4 (NP_073734.1).
Other Names	FRCP1
Accession No.	Swiss-Prot#:Q9H6D8NCBI Gene ID:64838
Uniprot	Q9H6D8
GeneID	64838;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

IHC 1:50 - 1:100


Images




Immunohistochemistry of paraffin-embedded human tonsil using FNDC4 at dilution of 1:150 (40x lens).




Immunohistochemistry of paraffin-embedded human colon using FNDC4 at dilution of 1:150 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using FNDC4 at dilution of 1:150 (40x lens).



Immunohistochemistry of paraffin-embedded human appendix using FNDC4 at dilution of 1:150 (40x lens).



Immunohistochemistry of paraffin-embedded mouse spleen using FNDC4 at dilution of 1:150 (40x lens).

Background

Note: This product is for in vitro research use only