

SPATS2 Polyclonal Antibody

Catalog No: #30166



Package Size: #30166-1 50ul #30166-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SPATS2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human SPATS2 (NP_075559.2).
Other Names	Nbla00526, P59SCR, SCR59, SPATA10
Accession No.	Swiss-Prot#:Q86XZ4NCBI Gene ID:65244
Uniprot	Q86XZ4
GeneID	65244;
Calculated MW	60kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Rat testis, using SPATS2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only