

PRKRIP1 Polyclonal Antibody

Catalog No: #30169



Package Size: #30169-1 50ul #30169-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PRKRIP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PRKRIP1 (NP_078929.1).
Other Names	C114, KRBOX3
Accession No.	Swiss-Prot#:Q9H875NCBI Gene ID:79706
Uniprot	Q9H875
GeneID	79706;
Calculated MW	21kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using PRKRIP1 at 1:1000 dilution.



Immunofluorescence analysis of C6 cells using PRKRIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of HeLa cells using PRKRIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of L929 cells using PRKRIP1 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only