

VPS37B Polyclonal Antibody

Catalog No: #30170

Package Size: #30170-1 50ul #30170-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	VPS37B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human VPS37B (NP_078943.1).
Accession No.	Swiss-Prot#:Q9H9H4NCBI Gene ID:79720
Uniprot	Q9H9H4
GeneID	79720;
Calculated MW	31kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Rat thymus, using VPS37B at 1:1000 dilution.

Background

Note: This product is for in vitro research use only