

ZNF768 Polyclonal Antibody

Catalog No: #30171



Package Size: #30171-1 50ul #30171-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF768 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human ZNF768 (NP_078947.3).
Accession No.	Swiss-Prot#:Q9H5H4NCBI Gene ID:79724
Uniprot	Q9H5H4
GeneID	79724;
Calculated MW	60kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse liver, using ZNF768 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only