

ST6GALNAC5 Polyclonal Antibody

Catalog No: #30176



Package Size: #30176-1 50ul #30176-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ST6GALNAC5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ST6GALNAC5 (NP_112227.1).
Other Names	SIAT7-E, SIAT7E, ST6GalNAcV
Accession No.	Swiss-Prot#:Q9BVH7NCBI Gene ID:81849
Uniprot	Q9BVH7
GeneID	81849;
Calculated MW	38kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ST6GALNAC5 at 1:1000 dilution.

Background

The protein encoded by this gene is a Golgi type II transmembrane glycosyltransferase. The encoded protein catalyzes the transfer of sialic acid to cell surface proteins to modulate cell-cell interactions. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016]

Note: This product is for in vitro research use only