

PUS7L Polyclonal Antibody

Catalog No: #30177



Package Size: #30177-1 50ul #30177-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PUS7L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human PUS7L (NP_112582.3).
Accession No.	Swiss-Prot#:Q9H0K6NCBI Gene ID:83448
Uniprot	Q9H0K6
GeneID	83448;
Calculated MW	81kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PUS7L at 1:1000 dilution.

Background

Note: This product is for in vitro research use only