

EAF1 Polyclonal Antibody

Catalog No: #30182



Package Size: #30182-1 50ul #30182-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	EAF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human EAF1 (NP_149074.3).
Accession No.	Swiss-Prot#:Q96JC9NCBI Gene ID:85403
Uniprot	Q96JC9
GeneID	85403;
Calculated MW	29kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using EAF1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only