

N4BP2L1 Polyclonal Antibody

Catalog No: #30185



Package Size: #30185-1 50ul #30185-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	N4BP2L1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human N4BP2L1 (NP_001073159).
Other Names	CG018
Accession No.	Swiss-Prot#:Q5TBK1NCBI Gene ID:90634
Uniprot	Q5TBK1
GeneID	90634;
Calculated MW	29kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse liver, using N4BP2L1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only