

DEPDC7 Polyclonal Antibody

Catalog No: #30187



Package Size: #30187-1 50ul #30187-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DEPDC7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human DEPDC7 (NP_001070710.1).
Other Names	TR2, dJ85M6.4
Accession No.	Swiss-Prot#:Q96QD5NCBI Gene ID:91614
Uniprot	Q96QD5
GeneID	91614;
Calculated MW	58kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DEPDC7 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only