

NAF1 Polyclonal Antibody

Catalog No: #30189



Package Size: #30189-1 50ul #30189-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NAF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NAF1 (NP_612395.2).
Accession No.	Swiss-Prot#:Q96HR8NCBI Gene ID:92345
Uniprot	Q96HR8
GeneID	92345;
Calculated MW	62-70kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NAF1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only