

OSBPL8 Polyclonal Antibody

Catalog No: #30194



Package Size: #30194-1 50ul #30194-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	OSBPL8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human OSBPL8 (NP_065892.1).
Other Names	MST120, MSTP120, ORP8, OSBP10
Accession No.	Swiss-Prot#:Q9BZF1NCBI Gene ID:114882
Uniprot	Q9BZF1
GeneID	114882;
Calculated MW	101kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using OSBPL8 at 1:1000 dilution.

Background

This gene encodes a member of a family of proteins containing an N-terminal pleckstrin homology domain and a highly conserved C-terminal oxysterol-binding protein-like sterol-binding domain. It binds multiple lipid-containing molecules, including phosphatidylserine, phosphatidylinositol 4-phosphate (PI4P) and oxysterol, and promotes their exchange between the endoplasmic reticulum and the plasma membrane. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]

Note: This product is for in vitro research use only