

SFXN2 Polyclonal Antibody

Catalog No: #30195



Package Size: #30195-1 50ul #30195-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SFXN2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SFXN2 (NP_849189.1).
Other Names	SFXN2
Accession No.	Uniprot:Q96NB2GeneID:118980
Uniprot	Q96NB2
GeneID	118980
Calculated MW	36kDa
SDS-PAGE MW	36KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using SFXN2 antibody.



Western blot analysis of extracts of various cell lines, using SFXN2 antibody.

Background

Mitochondrial amino-acid transporter that mediates transport of serine into mitochondria.

Note: This product is for in vitro research use only