

ASB17 Polyclonal Antibody

Catalog No: #30199

Package Size: #30199-1 50ul #30199-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ASB17 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ASB17 (NP_543144.1).
Other Names	Asb-17
Accession No.	Swiss-Prot#:Q8WXJ9NCBI Gene ID:127247
Uniprot	Q8WXJ9
GeneID	127247;
Calculated MW	34kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ASB17 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only