

NUDT16 Polyclonal Antibody

Catalog No: #30200



Package Size: #30200-1 50ul #30200-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NUDT16 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human NUDT16 (NP_689608.2).
Accession No.	Swiss-Prot#:Q96DE0NCBI Gene ID:131870
Uniprot	Q96DE0
GeneID	131870;
Calculated MW	21kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NUDT16 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only