

PNLDC1 Polyclonal Antibody

Catalog No: #30204



Package Size: #30204-1 50ul #30204-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PNLDC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant fusion protein of human PNLDC1 (NP_775787.1).
Other Names	HsPNLDC1
Accession No.	Swiss-Prot#:Q8NA58NCBI Gene ID:154197
Uniprot	Q8NA58
GeneID	154197;
Calculated MW	60kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PNLDC1 antibody.

Background

Note: This product is for in vitro research use only