

DNAJC18 Polyclonal Antibody

Catalog No: #30206



Package Size: #30206-1 50ul #30206-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	DNAJC18 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu,Rt
Immunogen Description	Recombinant fusion protein of human DNAJC18 (NP_689899.1).
Accession No.	Swiss-Prot#:Q9H819NCBI Gene ID:202052
Uniprot	Q9H819
GeneID	202052;
Calculated MW	42kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using DNAJC18 antibody.

Background

Note: This product is for in vitro research use only