

LIN9 Polyclonal Antibody

Catalog No: #30209

Package Size: #30209-1 50ul #30209-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LIN9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human LIN9 (NP_775106.2).
Other Names	BARA;BARPsv;Lin-9;TGS;TGS1;TGS2;LIN9
Accession No.	Uniprot:Q5TKA1GeneID:286826
Uniprot	Q5TKA1
GeneID	286826
Calculated MW	70kDa
SDS-PAGE MW	70kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of Mouse testis, using LIN9 antibody.

Background

This gene encodes a tumor suppressor protein that inhibits DNA synthesis and oncogenic transformation through association with the retinoblastoma 1 protein. The encoded protein also interacts with a complex of other cell cycle regulators to repress cell cycle-dependent gene expression in non-dividing cells. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul

Note: This product is for in vitro research use only