

VWC2 Polyclonal Antibody

Catalog No: #30210



Package Size: #30210-1 50ul #30210-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	VWC2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Ms,Rt
Immunogen Description	Recombinant fusion protein of Human VWC2.
Other Names	A930041G11Rik, G11, PSST739, UNQ739, cradin
Accession No.	Swiss-Prot#:_NCBI Gene ID:319922
GeneID	:319922
Calculated MW	35kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using VWC2 antibody.

Background

Note: This product is for in vitro research use only