

MTCO1 Polyclonal Antibody

Catalog No: #30224



Package Size: #30224-1 50ul #30224-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MTCO1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Hu,Ms,Rt
Immunogen Description	Recombinant protein of human MTCO1.
Other Names	COI; MTCO1; COX1
Accession No.	Swiss-Prot#:P00395NCBI Gene ID:4512
Uniprot	P00395
GeneID	4512;
Calculated MW	40kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using MTCO1 antibody.



Immunofluorescence analysis of NIH-3T3 cells using MTCO1 Polyclonal Antibody.

Background

Note: This product is for in vitro research use only