

NANOS3 Polyclonal Antibody

Catalog No: #30231



Package Size: #30231-1 50ul #30231-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NANOS3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Hu
Immunogen Description	Recombinant fusion protein of human NANOS3 (NP_001092092.1).
Other Names	NOS3; NANOS1L; ZC2HC12C
Accession No.	Swiss-Prot#:P60323NCBI Gene ID:342977
Uniprot	P60323
GeneID	342977;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only