

## NNAT Polyclonal Antibody

Catalog No: #30237



Package Size: #30237-1 50ul #30237-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                   |
|-----------------------|-------------------------------------------------------------------|
| Product Name          | NNAT Polyclonal Antibody                                          |
| Host Species          | Rabbit                                                            |
| Clonality             | Polyclonal                                                        |
| Isotype               | IgG                                                               |
| Purification          | Affinity purification                                             |
| Applications          | WB                                                                |
| Species Reactivity    | Hu                                                                |
| Immunogen Description | Recombinant protein of human NNAT.                                |
| Accession No.         | Swiss-Prot#:Q16517NCBI Gene ID:4826                               |
| Uniprot               | Q16517                                                            |
| GeneID                | 4826;                                                             |
| Calculated MW         | Refer to figures                                                  |
| Formulation           | Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4. |
| Storage               | Store at -20°C                                                    |

## Application Details

WB□1:500 - 1:2000

## Background

Note: This product is for in vitro research use only