

IRF2BP2 Polyclonal Antibody

Catalog No: #30245



Package Size: #30245-1 50ul #30245-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	IRF2BP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of human IRF2BP2.
Accession No.	Swiss-Prot#:Q7Z5L9NCBI Gene ID:359948
Uniprot	Q7Z5L9
GeneID	359948;
Calculated MW	61kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using IRF2BP2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only