

NUP214 Antibody

Catalog No: #46102

Package Size: #46102-1 50ul #46102-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NUP214 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	NUP214 Antibody detects endogenous levels of total NUP214
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NUP214
Target Name	NUP214
Other Names	NUP214, 214 kDa nucleoporin, CAN, CAN protein, putative oncogene, D9S46E, N214, Nucleoporin 214kD (CAIN), Nucleoporin 214kDa, Nucleoporin Nup214, p250, CAIN, KIAA0023, Protein CAN
Accession No.	Swiss-Prot#: P35658NCBI Gene ID:8021
Uniprot	P35658
GeneID	8021;
Calculated MW	214kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis NUP214 using LOVO whole cell lysates

Note: This product is for in vitro research use only