

UBE2S Polyclonal Antibody

Catalog No: #30283



Package Size: #30283-1 50ul #30283-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	UBE2S Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human UBE2S (NP_055316.2).
Other Names	UBE2S; E2-EPF; E2EPF; EPF5; ubiquitin-conjugating enzyme E2 S
Accession No.	Swiss-Prot#:Q16763NCBI Gene ID:27338
Uniprot	Q16763
GeneID	27338;
Calculated MW	26kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts from normal (control) and UBE2S knockout (KO) 293T cells, using UBE2S at 1:1000 dilution.

Background

This gene encodes a member of the ubiquitin-conjugating enzyme family. The encoded protein is able to form a thiol ester linkage with ubiquitin in a ubiquitin activating enzyme-dependent manner, a characteristic property of ubiquitin carrier proteins.

Note: This product is for in vitro research use only