

NLRX1 Polyclonal Antibody

Catalog No: #30284



Package Size: #30284-1 50ul #30284-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NLRX1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NLRX1 (NP_001269073.1).
Other Names	NLRX1; CLR11.3; DLNB26; NOD26; NOD5; NOD9; NLR family member X1
Accession No.	Swiss-Prot#:Q86UT6NCBI Gene ID:79671
Uniprot	Q86UT6
GeneID	79671;
Calculated MW	108kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts from normal (control) and NLRX1 knockout (KO) HeLa cells, using NLRX1 at 1:3000 dilution.

Background

The protein encoded by this gene is a member of the NLR family and localizes to the outer mitochondrial membrane. The encoded protein is a regulator of mitochondrial antiviral responses. Three transcript variants encoding the same protein have been found for this gene.

Note: This product is for in vitro research use only