

ZDHHC6 Polyclonal Antibody

Catalog No: #30310



Package Size: #30310-1 50ul #30310-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZDHHC6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ZDHHC6.
Accession No.	Swiss-Prot#:Q9H6R6NCBI Gene ID:64429
Uniprot	Q9H6R6
GeneID	64429;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Rat lung, using ZDHHC6 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only