

MEK1/MEK2 Polyclonal Antibody

Catalog No: #30313



Package Size: #30313-1 50ul #30313-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	MEK1/MEK2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human MAP2K2
Accession No.	Swiss-Prot#:Q02750/P36507NCBI Gene ID:5604/5605
Uniprot	Q02750
GeneID	5604;
Calculated MW	45kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of PC-12 cells, using MAP2K1/MAP2K2 at 1:500 dilution.



Western blot analysis of extracts of various cell lines, using MAP2K1/MAP2K2 at 1:500 dilution.

Background

Note: This product is for in vitro research use only