

FBXW8 Polyclonal Antibody

Catalog No: #30316

Package Size: #30316-1 50ul #30316-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FBXW8 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human FBXW8.
Accession No.	Swiss-Prot#:Q8N3Y1NCBI Gene ID:26259
Uniprot	Q8N3Y1
GeneID	26259;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:200

Images



Immunofluorescence analysis of L929 cells using FBXW8 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only