

SMG5 Polyclonal Antibody

Catalog No: #30317



Package Size: #30317-1 50ul #30317-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SMG5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant protein of human SMG5.
Accession No.	Swiss-Prot#:Q9UPR3NCBI Gene ID:23381
Uniprot	Q9UPR3
GeneID	23381;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only