

TBCCD1 Polyclonal Antibody

Catalog No: #30329



Package Size: #30329-1 50ul #30329-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TBCCD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human TBCCD1.
Other Names	TBCCD1
Accession No.	Uniprot:Q9NVR7GeneID:55171
Uniprot	Q9NVR7
GeneID	55171
Calculated MW	Refer to figures
SDS-PAGE MW	70kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB1:500 - 1:2000IHC1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using TBCCD1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Mouse kidney using TBCCD1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Human placenta using TBCCD1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded Rat kidney using TBCCD1 Rabbit pAb.

Background

Plays a role in the regulation of centrosome and Golgi apparatus positioning, with consequences on cell shape and cell migration.

Note: This product is for in vitro research use only