

K48-linkage Specific Polyubiquitin pAb

Catalog No: #30339



Package Size: #30339-1 50ul #30339-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	K48-linkage Specific Polyubiquitin pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide corresponding to a sequence containing polyubiquitinated protein (K48 linkage).
Calculated MW	Refer to figures
SDS-PAGE MW	10-100KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using K48-linkage Specific Polyubiquitin antibody.

Note: This product is for in vitro research use only