

FAM192A Polyclonal Antibody

Catalog No: #30344



Package Size: #30344-1 50ul #30344-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FAM192A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human, Mouse
Immunogen Description	Recombinant protein of human FAM192A.
Accession No.	Swiss-Prot#: Q9GZU8 NCBI Gene ID: 80011
Uniprot	Q9GZU8
GeneID	80011;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB □ 1:500 - 1:2000

Background

Note: This product is for in vitro research use only