

NOL4 Antibody

Catalog No: #46109



Package Size: #46109-1 50ul #46109-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NOL4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	NOL4 Antibody detects endogenous levels of total NOL4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NOL4
Target Name	NOL4
Other Names	NOL4, Cancer/testis antigen 125, CT125, HRIHFB2255, Nucleolar-localized protein, NOLP, Nucleolar localized protein, Nucleolar protein 4
Accession No.	Swiss-Prot#: O94818NCBI Gene ID:8715
Uniprot	O94818
GeneID	8715;
Calculated MW	71 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only