

SUPT3H Polyclonal Antibody

Catalog No: #30347



Package Size: #30347-1 50ul #30347-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SUPT3H Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human SUPT3H.
Accession No.	Swiss-Prot#:O75486NCBI Gene ID:8464
Uniprot	O75486
GeneID	8464;
Calculated MW	36kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of A-549 cells, using SUPT3H at 1:1000 dilution.

Background

Note: This product is for in vitro research use only