

Mrpl48 Polyclonal Antibody

Catalog No: #30350



Package Size: #30350-1 50ul #30350-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Mrpl48 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human Mrpl48 .
Other Names	L48mt; CGI-118; MRP-L48; D4Ert786e; 1810030E20Rik; 2610028L11Rik
Accession No.	Swiss-Prot#:Q8JZS9NCBI Gene ID:52443
Uniprot	Q8JZS9
GeneID	52443;
Calculated MW	19kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using Mrpl48 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only