

CFAP61 Polyclonal Antibody

Catalog No: #30353



Package Size: #30353-1 50ul #30353-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CFAP61 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human CFAP61.
Other Names	C20orf26;CaM-IP3;dJ1002M8.3;dJ1178H5.4;CFAP61
Accession No.	Uniprot:Q8NHU2GeneID:26074
Uniprot	Q8NHU2
GeneID	26074
Calculated MW	Refer to figures
SDS-PAGE MW	141KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images



Western blot analysis of extracts of Mouse eye, using CFAP61 antibody.



Immunohistochemistry of paraffin-embedded Human colon using CFAP61 Rabbit pAb.



Immunofluorescence analysis of Mouse testis using CFAP61 Rabbit pAb.

Background

May regulate cilium motility through its role in the assembly of the axonemal radial spokes.

Note: This product is for in vitro research use only