

K29-linkage Specific Polyubiquitin pAb

Catalog No: #30354



Package Size: #30354-1 50ul #30354-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	K29-linkage Specific Polyubiquitin pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human K29-linkage Specific Polyubiquitin.
Calculated MW	Refer to figures
SDS-PAGE MW	/
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using K29-linkage Specific Polyubiquitin Rabbit pAb.

Note: This product is for in vitro research use only