

Mrpl58 Polyclonal Antibody

Catalog No: #30361

Package Size: #30361-1 50ul #30361-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mrpl58 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human Mrpl58.
Other Names	Ict1; MRP-L58; 1110001A02Rik; 1110002E03Rik
Accession No.	Swiss-Prot#:Q8R035NCBI Gene ID:68572
Uniprot	Q8R035
GeneID	68572;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only