

IGF1R/INSR Polyclonal Antibody

Catalog No: #30365



Package Size: #30365-1 50ul #30365-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IGF1R/INSR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant protein of human IGF1R/INSR.
Accession No.	Swiss-Prot#:P08069/P06213NCBI Gene ID:3480/3643
Uniprot	P08069
GeneID	3480;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only