

RPGRIP1L Polyclonal Antibody

Catalog No: #30366



Package Size: #30366-1 50ul #30366-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RPGRIP1L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human RPGRIP1L.
Other Names	FTM; MKS5; CORS3; JBTS7; NPHP8; PPP1R134
Accession No.	Swiss-Prot#:Q68CZ1NCBI Gene ID:23322
Uniprot	Q68CZ1
GeneID	23322;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

The protein encoded by this gene can localize to the basal body-centrosome complex or to primary cilia and centrosomes in ciliated cells. The encoded protein has been found to interact with nephrocystin-4. Defects in this gene are a cause of Joubert syndrome type 7 (JBTS7) and Meckel syndrome type 5 (MKS5).

Note: This product is for in vitro research use only