

PHLPP2 Polyclonal Antibody

Catalog No: #30368



Package Size: #30368-1 50ul #30368-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PHLPP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human PHLPP2.
Other Names	PPM3B; PHLPPL
Accession No.	Swiss-Prot#:Q6ZVD8NCBI Gene ID:23035
Uniprot	Q6ZVD8
GeneID	23035;
Calculated MW	140kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using PHLPP2 at 1:1000 dilution.



Immunofluorescence analysis of L929 cells using PHLPP2 at dilution of 1:100. Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only