

RMND1 Polyclonal Antibody

Catalog No: #30380



Package Size: #30380-1 50ul #30380-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RMND1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human RMND1.
Other Names	RMD1; C6orf96; COXPD11; bA351K16; bA351K16.3
Accession No.	Swiss-Prot#:Q9NWS8NCBI Gene ID:55005
Uniprot	Q9NWS8
GeneID	55005;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

The protein encoded by this gene belongs to the evolutionary conserved sif2 family of proteins that share the DUF155 domain in common. This protein is thought to be localized in the mitochondria and involved in mitochondrial translation. Mutations in this gene are associated with combined oxidative phosphorylation deficiency-11. Alternatively spliced transcript variants have been found for this gene.

Note: This product is for in vitro research use only