

CHAF1A Polyclonal Antibody

Catalog No: #30412



Package Size: #30412-1 50ul #30412-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CHAF1A Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human CHAF1A (NP_005474.2).
Other Names	CHAF1A; CAF-1; CAF1; CAF1B; CAF1P150; P150; chromatin assembly factor 1 subunit A
Accession No.	Swiss-Prot#:Q13111NCBI Gene ID:10036
Uniprot	Q13111
GeneID	10036;
Calculated MW	160kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using CHAF1A at 1:1000 dilution.

Background

Note: This product is for in vitro research use only