

N5-Methylcytidine Polyclonal Antibody

Catalog No: #30415



Package Size: #30415-1 50ul #30415-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	N5-Methylcytidine Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human 5mC
Other Names	m5c
SDS-PAGE MW	/
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

IF 1:50 - 1:200

Images



Immunofluorescence analysis of HeLa cells using N5-Methylcytidine antibody.

Note: This product is for in vitro research use only