

STMN4 Polyclonal Antibody

Catalog No: #30430

Package Size: #30430-1 50ul #30430-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	STMN4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant fusion protein of human STMN4 (NP_110422.2).
Other Names	STMN4; RB3; stathmin-4
Accession No.	Swiss-Prot#:Q9H169NCBI Gene ID:81551
Uniprot	Q9H169
GeneID	81551;
Calculated MW	33kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:1000

Images



Western blot analysis of extracts of various cell lines, using STMN4 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only