

FAM117B Polyclonal Antibody

Catalog No: #30466



Package Size: #30466-1 50ul #30466-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM117B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human FAM117B (NP_775782.2).
Other Names	FAM117B; ALS2CR13; protein FAM117B
Accession No.	Swiss-Prot#:Q6P1L5NCBI Gene ID:150864
Uniprot	Q6P1L5
GeneID	150864;
Calculated MW	65kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using FAM117B .

Background

Note: This product is for in vitro research use only