

PPP3CB Polyclonal Antibody

Catalog No: #30473



Package Size: #30473-1 50ul #30473-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PPP3CB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human PPP3CB (NP_001135825.1).
Other Names	PPP3CB; CALNA2; CALNB; CNA2; PP2Bbeta; protein phosphatase 3 catalytic subunit beta
Accession No.	Swiss-Prot#:P16298NCBI Gene ID:5532
Uniprot	P16298
GeneID	5532;
Calculated MW	69kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse brain, using PPP3CB at 1:1000 dilution.

Background

Note: This product is for in vitro research use only