

ANO2 Polyclonal Antibody

Catalog No: #30478



Package Size: #30478-1 50ul #30478-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ANO2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ANO2 (NP_001265526.1).
Other Names	ANO2; C12orf3; TMEM16B; anoctamin-2
Accession No.	Swiss-Prot#:Q9NQ90NCBI Gene ID:57101
Uniprot	Q9NQ90
GeneID	57101;
Calculated MW	114kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:100

Images



Immunofluorescence analysis of Y79 cells using ANO2 Polyclonal at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Western blot analysis of extracts of mouse brain, using ANO2 at 1:3000 dilution.

Background

Note: This product is for in vitro research use only